



Anti-ZNF229 (aa 231-342) polyclonal antibody (DPABH-15239)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in transcriptional regulation.
Immunogen	Recombinant fragment corresponding to amino acids 231-342 of Human ZNF229 (Q9UJW7).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [ZNF229 zinc finger protein 230 \[Homo sapiens \]](#)

Official Symbol	ZNF229
Synonyms	ZNF229; zinc finger protein 229;
Entrez Gene ID	7772
Protein Refseq	NP_001265439.2
UniProt ID	A0A087WVZ7
Function	DNA binding; metal ion binding; sequence-specific DNA binding transcription factor activity;