



Anti-ZNF615 polyclonal antibody (DPABH-15249)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in transcriptional regulation.
Immunogen	Recombinant fragment from a region of Human ZNF615.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [ZNF615 zinc finger protein 616 \[Homo sapiens \]](#)

Official Symbol	ZNF615
Synonyms	ZNF615; zinc finger protein 615;
Entrez Gene ID	284370
Protein Refseq	NP_001186253.1
UniProt ID	Q8N8J6
Pathway	Gene Expression;
Function	DNA binding; metal ion binding; sequence-specific DNA binding transcription factor activity;