



Anti-RPGR (aa 379-509) polyclonal antibody (DPABH-15273)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Could be a guanine-nucleotide releasing factor.
Immunogen	antigen: corresponding to internal sequence amino acids 379-509 of the Human RPGR protein (Q92834).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [RPGR retinitis pigmentosa GTPase regulator \[Homo sapiens \]](#)

Official Symbol	RPGR
Synonyms	RPGR; retinitis pigmentosa GTPase regulator; CRD; RP3; COD1; PCDX; RP15; XLRP3; orf15; CORDX1; X-linked retinitis pigmentosa GTPase regulator; retinitis pigmentosa 15; retinitis pigmentosa 3 GTPase regulator;
Entrez Gene ID	6103
Protein Refseq	NP_000319.1
UniProt ID	Q92834
Function	guanyl-nucleotide exchange factor activity; poly(A) RNA binding; protein binding;