



# Anti-FAM181A (aa 232-352) polyclonal antibody (DPABH-15284)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	This gene encodes a protein that contains a coiled-coil domain. Several alternatively spliced transcript variants encoding multiple isoforms have been found for this gene.
<b>Immunogen</b>	antigen, corresponding to C-terminus amino acids 232-352 of the Human FAM181A protein (Q8N9Y4).
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ICC/IF, IHC-P
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">FAM181A family with sequence similarity 182, member A [ Homo sapiens ]</a>
<b>Official Symbol</b>	FAM181A
<b>Synonyms</b>	FAM181A; family with sequence similarity 181, member A; C14orf152; protein FAM181A;
<b>Entrez Gene ID</b>	<a href="#">90050</a>
<b>Protein Refseq</b>	<a href="#">NP_001194000.1</a>
<b>UniProt ID</b>	<a href="#">Q8N9Y4</a>