



Anti-KIAA2026 (aa 1844-1965) polyclonal antibody (DPABH-15315)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	KIAA2026 (KIAA2026) is a protein-coding gene.
Immunogen	Recombinant fragment: APGGTTIHTA SAPSNVTSLV GSQFSEPCIQ QKIVINTSTP LAPGTQIMIN GTRFIVPPQG LGAGSHVLLI STNPKYGAPL VLNSGQGIQS TPIDNSAQKI TLASNNSLSG QPLQHPLRSP TK, corresponding to C terminal amino acids 1844-1965 of Human KIAA2026 (UniProt ID: Q5HYC2).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	KIAA2026 KIAA2027 [Homo sapiens]
Official Symbol	KIAA2026
Synonyms	KIAA2026; uncharacterized protein KIAA2026;
Entrez Gene ID	158358
Protein Refseq	NP_001017969.2
UniProt ID	Q5HYC2