



# Anti-CXORF36 (aa 31-164) polyclonal antibody (DPABH-15336)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                            |                                                                                                                                                                                   |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Antigen Description</b> | CXorf36 (chromosome X open reading frame 36) is a protein-coding gene. Diseases associated with CXorf36 include mental retardation. An important paralog of this gene is C3orf58. |
| <b>Immunogen</b>           | Recombinant fragment corresponding to amino acids 31-164 of Human CXorf36 (Q9H7Y0).                                                                                               |
| <b>Isotype</b>             | IgG                                                                                                                                                                               |
| <b>Source/Host</b>         | Rabbit                                                                                                                                                                            |
| <b>Species Reactivity</b>  | Human                                                                                                                                                                             |
| <b>Purification</b>        | Immunogen affinity purified                                                                                                                                                       |
| <b>Conjugate</b>           | Unconjugated                                                                                                                                                                      |
| <b>Applications</b>        | IHC-P                                                                                                                                                                             |
| <b>Format</b>              | Liquid                                                                                                                                                                            |
| <b>Size</b>                | 100 µl                                                                                                                                                                            |
| <b>Buffer</b>              | pH: 7.20; Constituents: 59% PBS, 40% Glycerol                                                                                                                                     |
| <b>Preservative</b>        | 0.02% Sodium Azide                                                                                                                                                                |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.                                                                                    |

## GENE INFORMATION

|                        |                                                                                                                                                                                                                 |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>       | <a href="#">CXORF36 chromosome X open reading frame 37 [ Homo sapiens ]</a>                                                                                                                                     |
| <b>Official Symbol</b> | CXORF36                                                                                                                                                                                                         |
| <b>Synonyms</b>        | CXORF36; chromosome X open reading frame 36; DIA1R; PRO3743; EPQL1862; bA435K1.1; 4930578C19Rik; deleted in autism-related protein 1; hCG1981635; UPF0672 protein CXorf36; deleted in autism-1 related protein; |
| <b>Entrez Gene ID</b>  | <a href="#">79742</a>                                                                                                                                                                                           |
| <b>Protein Refseq</b>  | <a href="#">NP_078965.2</a>                                                                                                                                                                                     |
| <b>UniProt ID</b>      | <a href="#">Q9H7Y0</a>                                                                                                                                                                                          |