



Anti-C6ORF211 (aa 52-172) polyclonal antibody (DPABH-15396)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C6orf211 (chromosome 6 open reading frame 211) is a protein-coding gene. Diseases associated with C6orf211 include breast cancer susceptibility , and cleft lip.
Immunogen	Recombinant fragment, corresponding to amino acids 52-172 of Human C6orf211.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Store at -20°C.

GENE INFORMATION

Gene Name	C6ORF211 chromosome 6 open reading frame 212 [Homo sapiens]
------------------	---

Official Symbol	C6ORF211
Synonyms	C6ORF211; chromosome 6 open reading frame 211; UPF0364 protein C6orf211;
Entrez Gene ID	79624
Protein Refseq	NP_001273491.1
UniProt ID	B4DPT6