



Anti-LHFPL1 (aa 39-85) polyclonal antibody (DPABH-15405)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	LHFPL1 (lipoma HMGIC fusion partner-like 1) is a protein-coding gene. Diseases associated with LHFPL1 include lipoma. An important paralog of this gene is LHFPL5.
Immunogen	antigen sequence, corresponding to amino acids 39-85 (GKPVSFSTFR RCNYPVRGEG HSLIMVEECG RYASFNAIPS LAWQMCT) of Human LHFPL1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	LHFPL1 lipoma HMGIC fusion partner-like 2 [Homo sapiens]
Official Symbol	LHFPL1
Synonyms	LHFPL1; lipoma HMGIC fusion partner-like 1; lipoma HMGIC fusion partner-like 1 protein; LHP-like protein 1;
Entrez Gene ID	340596
Protein Refseq	NP_835469.1
UniProt ID	Q86WI0