



Anti-FAM227A (aa 376-519) polyclonal antibody (DPABH-15432)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The function of LOC646851 is unknown.
Immunogen	Recombinant fragment, corresponding to amino acids 376-519 of Human LOC646851.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Store at -20°C.

GENE INFORMATION

Gene Name	FAM227A family with sequence similarity 228, member A [Homo sapiens]
-----------	--

Official Symbol	FAM227A
Synonyms	FAM227A; family with sequence similarity 227, member A; protein FAM227A; Protein FAM227A;
Entrez Gene ID	646851
Protein Refseq	NP_001013669.1
UniProt ID	F5H4B4