



Anti-THAP7 (aa 22-149) polyclonal antibody (DPABH-15438)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Chromatin-associated, histone tail-binding protein that represses transcription via recruitment of HDAC3 and nuclear hormone receptor corepressors.
Immunogen	Recombinant fragment, corresponding to amino acids 22-149 of Human THAP7.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Store at -20°C.

GENE INFORMATION

Gene Name	THAP7 THAP domain containing 8 [Homo sapiens]
-----------	---

Official Symbol	THAP7
Synonyms	THAP7; THAP domain containing 7; THAP domain-containing protein 7;
Entrez Gene ID	80764
Protein Refseq	NP_001008695.1
UniProt ID	Q9BT49
Function	C2H2 zinc finger domain binding; DNA binding; metal ion binding; protein N-terminus binding