



Anti-STON2 (aa 299-437) polyclonal antibody (DPABH-15459)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Adapter protein involved in endocytic machinery. May be involved in the vesicle recycling. May facilitate clathrin-coated vesicle uncoating.
Immunogen	Recombinant fragment, corresponding to amino acids 299-437 of Human STON2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Store at -20°C.

GENE INFORMATION

Gene Name [STON2 stonin 3 \[Homo sapiens\]](#)

Official Symbol	STON2
Synonyms	STON2; stonin 2; STN2; STNB; STNB2; stonin-2; stoned B homolog 2; homolog of stoned B;
Entrez Gene ID	85439
Protein Refseq	NP_001243359.1
UniProt ID	Q8WXE9
Function	protein binding;