



# Anti-CCDC176 (aa 1-95) polyclonal antibody (DPABH-15470)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                     |                                                                            |
|---------------------|----------------------------------------------------------------------------|
| Antigen Description | CCDC176 (coiled-coil domain containing 176) is a protein-coding gene.      |
| Immunogen           | Recombinant fragment, corresponding to amino acids 1-95 of Human C14orf45. |
| Isotype             | IgG                                                                        |
| Source/Host         | Rabbit                                                                     |
| Species Reactivity  | Human                                                                      |
| Purification        | Immunogen affinity purified                                                |
| Conjugate           | Unconjugated                                                               |
| Applications        | IHC-P                                                                      |
| Format              | Liquid                                                                     |
| Size                | 100 µl                                                                     |
| Buffer              | pH: 7.20; Constituents: 59% PBS, 40% Glycerol                              |
| Preservative        | 0.02% Sodium Azide                                                         |
| Storage             | Shipped at 4°C. Store at -20°C.                                            |

## GENE INFORMATION

|           |                                                                            |
|-----------|----------------------------------------------------------------------------|
| Gene Name | <a href="#">CCDC176 coiled-coil domain containing 177 [ Homo sapiens ]</a> |
|-----------|----------------------------------------------------------------------------|

|                        |                                                                                                                                   |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| <b>Official Symbol</b> | CCDC176                                                                                                                           |
| <b>Synonyms</b>        | CCDC176; coiled-coil domain containing 176; C14orf45; basal body-orientation factor 1; coiled-coil domain-containing protein 176; |
| <b>Entrez Gene ID</b>  | <a href="#">80127</a>                                                                                                             |
| <b>Protein Refseq</b>  | <a href="#">NP_079333.2</a>                                                                                                       |
| <b>UniProt ID</b>      | <a href="#">Q8ND07</a>                                                                                                            |