



Anti-CCDC176 (aa 376-522) polyclonal antibody (DPABH-15507)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CCDC176 (coiled-coil domain containing 176) is a protein-coding gene.
Immunogen	Recombinant fragment, corresponding to amino acids 376-522 of Human C14orf45 (Q8ND07).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [CCDC176 coiled-coil domain containing 177 \[Homo sapiens \]](#)

Official Symbol	CCDC176
Synonyms	CCDC176; coiled-coil domain containing 176; C14orf45; basal body-orientation factor 1; coiled-coil domain-containing protein 176;
Entrez Gene ID	80127
Protein Refseq	NP_079333.2
UniProt ID	Q8ND07