



Anti-ZNF184 (aa 71-212) polyclonal antibody (DPABH-15552)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in transcriptional regulation.
Immunogen	Recombinant fragment, corresponding to amino acids 71-212 of Human ZNF184
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [ZNF184 zinc finger protein 185 \[Homo sapiens \]](#)

Official Symbol	ZNF184
Synonyms	ZNF184; zinc finger protein 184; zinc finger protein 184 (Kruppel-like);
Entrez Gene ID	7738
Protein Refseq	NP_009080.2
UniProt ID	A0A024RCM2
Pathway	Gene Expression.
Function	DNA binding; zinc ion binding;