



Anti-HYPM (aa 1-96) polyclonal antibody (DPABH-15617)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a WW-domain containing protein similar to yeast splicing factor PRP40. This protein has been shown to interact with Huntingtin and MeCP2. Alternative splicing results in different transcriptional variants.
Immunogen	antigen sequence: NPANIMSEKK NCKNSSTNNN QTQDPSRNEL QVPRSFVDRV VQDERDVQSQ SSSTINTLLT LLDCLADYIM ERVGLEASNN GSMRNTSQDR EREVDNNREP H, corresponding to amino acids 1-96 of Human CXorf27
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	HYPM huntingtin interacting protein M [Homo sapiens]
Official Symbol	HYPM
Synonyms	HYPM; huntingtin interacting protein M; HIP17; CXorf27; huntingtin-interacting protein M; huntingtin yeast partner M;
Entrez Gene ID	25763
Protein Refseq	NP_036406.1
UniProt ID	O75409
Function	protein heterodimerization activity;