



Anti-IFT43 (aa 107-202) polyclonal antibody (DPABH-15648)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Component of IFT complex A (IFT-A) involved in retrograde ciliary transport along microtubules from the ciliary tip to the base.
Immunogen	Recombinant fragment, corresponding to amino acids 107-202 (EEVQEEDFVL QVAAPPSIQI KRVMTYRDLD NDLMKYSAIQ TLDGEIDLKL LTKVLAPEHE VREDDVGWDW DHLFTEVSSE VLTEWDPLQT EKEDPA) of Human IFT43 (Q96FT9).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	IFT43 intraflagellar transport 44 homolog (Chlamydomonas) [Homo sapiens]
Official Symbol	IFT43
Synonyms	IFT43; intraflagellar transport 43 homolog (Chlamydomonas); CED3; C14orf179; intraflagellar transport protein 43 homolog; IFT complex A subunit;
Entrez Gene ID	112752
Protein Refseq	NP_001096034.1
UniProt ID	Q96FT9