



Anti-TMEM260 (aa 239-319) polyclonal antibody (DPABH-15664)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMEM260 (transmembrane protein 260) is a protein-coding gene.
Immunogen	antigen sequence corresponding amino acids 239-319 (SSYLNHARWT WGDQTTLQGF LTHFLREEYG TFSLAKSEIG SSMSEILLSQ VTNMRTELSF NIQALAVCAN ICLATKDRQN P) of Human C14orf101.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	TMEM260 transmembrane protein 261 [Homo sapiens]
Official Symbol	TMEM260
Synonyms	TMEM260; transmembrane protein 260; C14orf101; UPF0679 protein C14orf101;
Entrez Gene ID	54916
Protein Refseq	NP_060269.3
UniProt ID	B3KN73