



Anti-ASB1 (aa 106-244) polyclonal antibody (DPABH-15670)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May play a role in testis development (By similarity). Probable substrate-recognition component of a SCF-like ECS (Elongin-Cullin-SOCS-box protein) E3 ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.
Immunogen	Recombinant fragment, corresponding to amino acids 106-244 of Human ASB1 (Q9Y576).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	ASB1 ankyrin repeat and SOCS box containing 2 [Homo sapiens]
Official Symbol	ASB1
Synonyms	ASB1; ankyrin repeat and SOCS box containing 1; ASB-1; ankyrin repeat and SOCS box protein 1; ankyrin repeat and SOCS box-containing 1;
Entrez Gene ID	51665
Protein Refseq	NP_001035535.1
UniProt ID	Q9Y576
Pathway	Adaptive Immune System; Class I MHC mediated antigen processing & presentation;
Function	molecular_function;