



Anti-FSCB (aa 721-807) polyclonal antibody (DPABH-15681)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in the later stages of fibrous sheath biogenesis. Binds calcium.
Immunogen	Recombinant fragment, corresponding to amino acids 721-807 (ADLLLTEEFP IGEASAEVSP PPSEQTPED ALVENVSTEF QSPQVAGIPA VKLGSVVLEG EAKFEEVSKI NSVLKDLSNT NDGQAPT) of Human FSCB (Q5H9T9).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	FSCB fibrous sheath CA [Homo sapiens]
Official Symbol	FSCB
Synonyms	FSCB; fibrous sheath CA; C14orf155; fibrous sheath CABYR-binding protein;
Entrez Gene ID	84075
Protein Refseq	NP_115511.3
UniProt ID	Q5H9T9
Function	calcium ion binding;