



Anti-FCGBP (aa 289-417) polyclonal antibody (DPABH-15692)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in the maintenance of the mucosal structure as a gel-like component of the mucosa.
Immunogen	antigen sequence, corresponding to internal sequence amino acids 289-417 of Human FCGBP.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [FCGBP Fc fragment of IgG binding protein \[Homo sapiens \]](#)

Official Symbol	FCGBP
Synonyms	FCGBP; Fc fragment of IgG binding protein; FC(GAMMA)BP; IgGFc-binding protein; fcgammaBP; Human Fc gamma BP; IgG Fc binding protein; fcgamma-binding protein antigen;
Entrez Gene ID	8857
Protein Refseq	NP_003881.2
UniProt ID	Q9Y6R7
Function	protein binding;