



Anti-ZDHHC15 (aa 238-322) polyclonal antibody (DPABH-15701)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Palmitoyltransferase specific for GAP43 and DLG4/PSD95.
Immunogen	antigen sequence: VSRNKTLEA FCTPVFTSGP EKNGFNLGFI KNIQQVFGDK KKFWLPIGS SPGDGHSFPM RSMNESQNPL LANEETWEDN EDDNQ, corresponding to internal sequence amino acids 238-322 of Human ZDHHC15.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	ZDHHC15 zinc finger, DHHC-type containing 16 [Homo sapiens]
Official Symbol	ZDHHC15
Synonyms	ZDHHC15; zinc finger, DHHC-type containing 15; MRX91; palmitoyltransferase ZDHHC15; DHHC-15; zinc finger, DHHC domain containing 15; zinc finger DHHC domain-containing protein 15;
Entrez Gene ID	158866
Protein Refseq	NP_001139728.1
UniProt ID	B3KY34
Function	palmitoyltransferase activity; protein-cysteine S-palmitoyltransferase activity; zinc ion binding;