



Anti-C14ORF119 (aa 15-131) polyclonal antibody (DPABH-15713)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C14orf119 (chromosome 14 open reading frame 119) is a protein-coding gene.
Immunogen	antigen sequence, corresponding to amino acids 15-131 of Human C14orf119.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	C14ORF119 chromosome 14 open reading frame 120 [Homo sapiens]
-----------	---

Official Symbol	C14ORF119
Synonyms	C14ORF119; chromosome 14 open reading frame 119; uncharacterized protein C14orf119; My028 protein;
Entrez Gene ID	55017
Protein Refseq	NP_060394.1
UniProt ID	Q9NWQ9