



# Anti-ERICH1 (aa 300-435) polyclonal antibody (DPABH-15721)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                            |                                                                                                                                                    |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Antigen Description</b> | ERICH1 (glutamate-rich 1) is a protein-coding gene. Diseases associated with ERICH1 include pancreatic ductal adenocarcinoma , and adenocarcinoma. |
| <b>Immunogen</b>           | antigen sequence, corresponding to C terminal amino acids 300-435 of Human ERICH1.                                                                 |
| <b>Isotype</b>             | IgG                                                                                                                                                |
| <b>Source/Host</b>         | Rabbit                                                                                                                                             |
| <b>Species Reactivity</b>  | Human                                                                                                                                              |
| <b>Purification</b>        | Immunogen affinity purified                                                                                                                        |
| <b>Conjugate</b>           | Unconjugated                                                                                                                                       |
| <b>Applications</b>        | IHC-P, WB, ICC/IF                                                                                                                                  |
| <b>Format</b>              | Liquid                                                                                                                                             |
| <b>Size</b>                | 100 µl                                                                                                                                             |
| <b>Buffer</b>              | pH: 7.20; Constituents: 59% PBS, 40% Glycerol                                                                                                      |
| <b>Preservative</b>        | 0.02% Sodium Azide                                                                                                                                 |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.                                                     |

## GENE INFORMATION

Gene Name [ERICH1 glutamate-rich 2 \[ Homo sapiens \]](#)

|                        |                                                              |
|------------------------|--------------------------------------------------------------|
| <b>Official Symbol</b> | ERICH1                                                       |
| <b>Synonyms</b>        | ERICH1; glutamate-rich 1; HSPC319; glutamate-rich protein 1; |
| <b>Entrez Gene ID</b>  | <a href="#">157697</a>                                       |
| <b>Protein Refseq</b>  | <a href="#">NP_997215.1</a>                                  |
| <b>UniProt ID</b>      | <a href="#">Q86X53</a>                                       |