



Anti-AP1G2 (aa 586-708) polyclonal antibody (DPABH-15773)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May function in protein sorting in late endosomes or multivesicular bodies (MVBs). Involved in MVB-assisted maturation of hepatitis B virus (HBV).
Immunogen	antigen sequence, corresponding to amino acids 586-708 of Human AP1G2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [AP1G2 adaptor-related protein complex 1, gamma 3 subunit \[Homo sapiens \]](#)

Official Symbol	AP1G2
Synonyms	AP1G2; adaptor-related protein complex 1, gamma 2 subunit; G2AD; AP-1 complex subunit gamma-like 2; gamma2-adaptin; clathrin-associated/assembly/adaptor protein, large, gamma-2;
Entrez Gene ID	8906
Protein Refseq	NP_001269403.1
UniProt ID	Q86U03
Pathway	Lysosome;
Function	protein binding; protein transporter activity;