



Anti-TPPP2 (aa 8-150) polyclonal antibody (DPABH-15793)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May bind tubulin but has no microtubule bundling activity.
Immunogen	antigen sequence, corresponding to amino acids 8-150 of Human TPPP2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	TPPP2 tubulin polymerization-promoting protein family member 3 [Homo sapiens]
-----------	---

Official Symbol	TPPP2
Synonyms	TPPP2; tubulin polymerization-promoting protein family member 2; P18; CT152; C14orf8; p25beta; TPPP/p18; protein p25-beta;
Entrez Gene ID	122664
Protein Refseq	NP_776245.2
UniProt ID	P59282