



Anti-C1ORF198 (aa 164-309) polyclonal antibody (DPABH-15840)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C1orf198 (chromosome 1 open reading frame 198) is a protein-coding gene.
Immunogen	antigen sequence, corresponding to amino acids 164-309 of Human C1orf198.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, ICC/IF, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	C1ORF198 chromosome 1 open reading frame 199 [Homo sapiens]
-----------	---

Official Symbol	C1ORF198
Synonyms	C1ORF198; chromosome 1 open reading frame 198; uncharacterized protein C1orf198;
Entrez Gene ID	84886
Protein Refseq	NP_001129966.1
UniProt ID	Q9H425