



Anti-KIAA1217 (aa 794-898) polyclonal antibody (DPABH-15862)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Required for normal development of intervertebral disks.
Immunogen	antigen, corresponding to amino acids 794-898 of Human SKT (UniProt ID: Q5T5P2).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [KIAA1217 KIAA1218 \[Homo sapiens \]](#)

Official Symbol	KIAA1217
Synonyms	KIAA1217; SKT; sickle tail protein homolog; likely orthologue of Mus musculus enhancer trap locus 4;
Entrez Gene ID	56243
Protein Refseq	NP_001091970.1
UniProt ID	Q5T5P2
Function	molecular_function;