



Anti-LANCL3 (aa 72-157) polyclonal antibody (DPABH-15874)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	LANCL3 (LanC lantibiotic synthetase component C-like 3 (bacterial)) is a protein-coding gene. GO annotations related to this gene include catalytic activity. An important paralog of this gene is LANCL2.
Immunogen	Recombinant fragment: VAYMLYHVSQ SPLFATARER YLRS AKRLID ACARAE EWGE PDADTRAAFL LGGAGVYAVA TLVYHALGRS DYVQPLGKFR ALCAVC, corresponding to amino acids 72-157 of Human LANCL3 (UniProt Q6ZV70).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Store at -20°C.

GENE INFORMATION

Gene Name	LANCL3 LanC lantibiotic synthetase component C-like 4 (bacterial) [Homo sapiens]
Official Symbol	LANCL3
Synonyms	LANCL3; LanC lantibiotic synthetase component C-like 3 (bacterial); lanC-like protein 3;
Entrez Gene ID	347404
Protein Refseq	NP_001163802.1
UniProt ID	Q6ZV70
Function	catalytic activity;