



Rabbit anti-Human CHADL polyclonal antibody (DPABH-15902)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	antigen sequence, corresponding to amino acids 493-620 of Human CHADL.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	CHADL chondroadherin-like [Homo sapiens]
Official Symbol	CHADL

Synonyms	CHADL; chondroadherin-like; SLRR4B; chondroadherin-like protein;
Entrez Gene ID	150356
Protein Refseq	NP_612490.1
UniProt ID	Q6NUI6