



Anti-TCP11L1 (aa 2-140) polyclonal antibody (DPABH-15932)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene is the central gene in a cluster of three vanin genes on chromosome 6q23-q24. Extensive alternative splicing has been described; the two most common variants are represented as RefSeqs.
Immunogen	Recombinant fragment corresponding to amino acids 2-140 of Human TCP11L1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Store at -20°C.

GENE INFORMATION

Gene Name	TCP11L1 t-complex 11, testis-specific-like 2 [Homo sapiens]
Official Symbol	TCP11L1
Synonyms	TCP11L1; t-complex 11, testis-specific-like 1; dJ85M6.3; T-complex protein 11-like protein 1; t-complex 11-like 1;
Entrez Gene ID	55346
Protein Refseq	NP_001139013.1
UniProt ID	Q9NUJ3