



Anti-USP40 (aa 755-900) polyclonal antibody (DPABH-15959)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be catalytically inactive.
Immunogen	Recombinant fragment, corresponding to amino acids 755-900 of Human USP40.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Store at -20°C.

GENE INFORMATION

Gene Name [USP40 ubiquitin specific peptidase 41 \[Homo sapiens \]](#)

Official Symbol	USP40
Synonyms	USP40; ubiquitin specific peptidase 40; ubiquitin carboxyl-terminal hydrolase 40; ubiquitin thioesterase 40; deubiquitinating enzyme 40; ubiquitin thiolesterase 40; ubiquitin specific protease 40; ubiquitin-specific-processing protease 40;
Entrez Gene ID	55230
Protein Refseq	NP_060688.1
UniProt ID	Q9NVE5
Function	cysteine-type peptidase activity; ubiquitinyl hydrolase activity;