



Anti-CCDC88C (aa 1571-1698) polyclonal antibody (DPABH-15967)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Negative regulator of the canonical Wnt signaling pathway, acting downstream of DVL to inhibit CTNNB1/Beta-catenin stabilization.
Immunogen	Recombinant fragment corresponding to amino acids 1571-1698 of Human CCDC88C.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [CCDC88C coiled-coil domain containing 89C \[Homo sapiens \]](#)

Official Symbol	CCDC88C
Synonyms	CCDC88C; coiled-coil domain containing 88C; DAPLE; HKRP2; KIAA1509; protein Daple; hook-related protein 2; Dvl-associating protein with a high frequency of leucine residues;
Entrez Gene ID	440193
Protein Refseq	NP_001073883.2
UniProt ID	B4DZB8
Function	PDZ domain binding; protein self-association;