



# Anti-TEX33 (aa 101-233) polyclonal antibody (DPABH-16011)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                     |                                                                                       |
|---------------------|---------------------------------------------------------------------------------------|
| Antigen Description | TEX33 (testis expressed 33) is a protein-coding gene.                                 |
| Immunogen           | Recombinant fragment corresponding to amino acids 101-233 of Human EAN57 (O43247)     |
| Isotype             | IgG                                                                                   |
| Source/Host         | Rabbit                                                                                |
| Species Reactivity  | Human                                                                                 |
| Purification        | Immunogen affinity purified                                                           |
| Conjugate           | Unconjugated                                                                          |
| Applications        | ICC/IF, IHC-P, WB                                                                     |
| Format              | Liquid                                                                                |
| Size                | 100 µl                                                                                |
| Buffer              | pH: 7.20; Constituents: 59% PBS, 40% Glycerol                                         |
| Preservative        | 0.02% Sodium Azide                                                                    |
| Storage             | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

## GENE INFORMATION

Gene Name [TEX33 testis expressed 34 \[ Homo sapiens \]](#)

|                        |                                                                                                             |
|------------------------|-------------------------------------------------------------------------------------------------------------|
| <b>Official Symbol</b> | TEX33                                                                                                       |
| <b>Synonyms</b>        | TEX33; testis expressed 33; EAN57; C22orf33; cE81G9.2; testis-expressed sequence 33 protein; protein EAN57; |
| <b>Entrez Gene ID</b>  | <a href="#">339669</a>                                                                                      |
| <b>Protein Refseq</b>  | <a href="#">NP_001157329.1</a>                                                                              |
| <b>UniProt ID</b>      | <a href="#">O43247</a>                                                                                      |