



Anti-EGFL6 (aa 402-546) polyclonal antibody (DPABH-16017)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May bind integrin alpha-8/beta-1 and play a role in hair follicle morphogenesis. Promotes matrix assembly.
Immunogen	Recombinant fragment corresponding to amino acids 402-546 of Human EGFL6
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [EGFL6 EGF-like-domain, multiple 7 \[Homo sapiens \]](#)

Official Symbol	EGFL6
Synonyms	EGFL6; EGF-like-domain, multiple 6; W80; MAEG; epidermal growth factor-like protein 6; EGF-like protein 6; MAM and EGF domain containing; EGF repeat-containing protein 6; MAM and EGF domains-containing gene protein;
Entrez Gene ID	25975
Protein Refseq	NP_001161362.1
UniProt ID	Q8IUX8
Function	calcium ion binding; integrin binding;