



Anti-SCIMP (aa 47-141) polyclonal antibody (DPABH-16049)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	SCIMP (SLP adaptor and CSK interacting membrane protein) is a protein-coding gene. Diseases associated with SCIMP include fainting. GO annotations related to this gene include protein binding.
Immunogen	Recombinant fragment: RRGKKWEIAK PLKHKQVDEE KMYENVLNES PVQLPPLPPR NWPSLEDSSP QEAPSQPPAT YSLVNKVKNK KTVSIPSYIE PEDDYDDVEI PANTE, corresponding to amino acids 47-141 of Human C17orf87
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	SCIMP SLP adaptor and CSK interacting membrane protein [Homo sapiens]
Official Symbol	SCIMP
Synonyms	SCIMP; SLP adaptor and CSK interacting membrane protein; UNQ5783; C17orf87; SLP adapter and CSK-interacting membrane protein; DTFT5783; transmembrane protein C17orf87; SLP65/SLP76, Csk-interacting membrane protein;
Entrez Gene ID	388325
Protein Refseq	NP_001258771.1
UniProt ID	Q6UWF3
Function	protein binding;