



Anti-ISY1 (aa 23-146) polyclonal antibody (DPABH-16090)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May play a role in pre-mRNA splicing.
Immunogen	Recombinant fragment corresponding to amino acids 23-146 of Human ISY1 (Q9ULR0).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [ISY1 ISY2 splicing factor homolog \(S. cerevisiae\) \[Homo sapiens \]](#)

Official Symbol	ISY1
Synonyms	ISY1; ISY1 splicing factor homolog (S. cerevisiae); FSAP33; pre-mRNA-splicing factor ISY1 homolog; functional spliceosome-associated protein 33;
Entrez Gene ID	57461
Protein Refseq	NP_001186398.1
UniProt ID	Q9ULR0
Pathway	Spliceosome; Spliceosome, 35S U5-snRNP;
Function	poly(A) RNA binding; protein binding;