



Anti-ZSCAN22 (aa 128-139) polyclonal antibody (DPABH-16110)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in transcriptional regulation.
Immunogen	Recombinant fragment corresponding to amino acids 128-139 of Human ZNF50 (P10073).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [ZSCAN22 zinc finger and SCAN domain containing 23 \[Homo sapiens \]](#)

Official Symbol	ZSCAN22
Synonyms	ZSCAN22; zinc finger and SCAN domain containing 22; HKR2; ZNF50; zinc finger and SCAN domain-containing protein 22; zinc finger protein 50; GLI-Kruppel family member HKR2; krueppel-related zinc finger protein 2;
Entrez Gene ID	342945
Protein Refseq	NP_862829.1
UniProt ID	P10073
Function	DNA binding; metal ion binding; molecular_function; sequence-specific DNA binding transcription factor activity