



Anti-USP42 (aa 454-580) polyclonal antibody (DPABH-16131)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Deubiquitinating enzyme which may play an important role during spermatogenesis.
Immunogen	Recombinant fragment corresponding to amino acids 454-580 of Human USP42 (Q9H9J4).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [USP42 ubiquitin specific peptidase 43 \[Homo sapiens \]](#)

Official Symbol	USP42
Synonyms	USP42; ubiquitin specific peptidase 42; ubiquitin carboxyl-terminal hydrolase 42; ubiquitin thioesterase 42; deubiquitinating enzyme 42; ubiquitin thiolesterase 42; ubiquitin specific protease 42; ubiquitin-specific-processing protease 42;
Entrez Gene ID	84132
Protein Refseq	NP_115548.1
UniProt ID	Q9H9J4
Function	protein binding; ubiquitin-specific protease activity;