



Anti-CNST (aa 302-446) polyclonal antibody (DPABH-16149)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Required for targeting of connexins to the plasma membrane.
Immunogen	antigen, corresponding to amino acids 302-446 of Human Consortin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [CNST consortin, connexin sorting protein \[Homo sapiens \]](#)

Official Symbol	CNST
Synonyms	CNST; consortin, connexin sorting protein; C1orf71; PPP1R64; consortin; protein phosphatase 1, regulatory subunit 64;
Entrez Gene ID	163882
Protein Refseq	NP_001132931.1
UniProt ID	Q6PJW8
Function	connexin binding; phosphatase binding; protein binding;