



Rabbit anti-Human AFAP1L1 polyclonal antibody (DPABH-16166)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	AFAP1L1 (actin filament associated protein 1-like 1) is a protein-coding gene. Diseases associated with AFAP1L1 include spindle cell sarcoma , and sarcoma. GO annotations related to this gene include phospholipid binding. An important paralog of this gene is AFAP1L2.
Immunogen	antigen, corresponding to amino acids 525-659 of Human AFAP1L1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	AFAP1L1 actin filament associated protein 1-like 2 [Homo sapiens]
Official Symbol	AFAP1L1
Synonyms	AFAP1L1; actin filament associated protein 1-like 1; actin filament-associated protein 1-like 1; AFAP1-like protein 1;
Entrez Gene ID	134265
Protein Refseq	NP_001139809.1
UniProt ID	Q8TED9