



Anti-SELK (aa 48-84) polyclonal antibody (DPABH-16312)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	antigen sequence: KRRSYGNSSD SRYDDGRGPP GNPPRRMGRI NHLRGPS, corresponding to amino acids 48-84 of Human SELK.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	SELK is a protein-coding gene. Diseases associated with SELK include staphyloenterotoxemia , and mastitis.
---------------------	--

Keywords SELK; selenoprotein K; SelK; HSPC030; HSPC297;

GENE INFORMATION

Entrez Gene ID [58515](#)

UniProt ID [Q9Y6D0](#)
