



Anti-DPY19L3 (aa 557-691) polyclonal antibody (DPABH-16346)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	DPY19L3 (dpy-19-like 3 (C. elegans)) is a protein-coding gene. Diseases associated with DPY19L3 include bipolar disorder. An important paralog of this gene is DPY19L4.
Immunogen	antigen, corresponding to amino acids 557-691 of Human DPY19L3 (Q6ZPD9).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [DPY19L3 dpy-19-like 4 \(C. elegans\) \[Homo sapiens \]](#)

Official Symbol	DPY19L3
Synonyms	DPY19L3; dpy-19-like 3 (C. elegans); probable C-mannosyltransferase DPY19L3; dpy-19-like protein 3; protein dpy-19 homolog 3;
Entrez Gene ID	147991
Protein Refseq	NP_001166245.1
UniProt ID	Q6ZPD9
Function	transferase activity, transferring glycosyl groups;