



Anti-KIAA0195 (aa 429-556) polyclonal antibody (DPABH-16357)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	KIAA0195 (KIAA0195) is a protein-coding gene.
Immunogen	antigen, corresponding to amino acids 429-556 of Human KIAA0195 (Q12767).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	KIAA0195 KIAA0196 [Homo sapiens]
-----------	--

Official Symbol	KIAA0195
Synonyms	KIAA0195; TMEM94; uncharacterized protein KIAA0195; transmembrane protein 94;
Entrez Gene ID	9772
Protein Refseq	NP_055553.3
UniProt ID	A0A024R8S6