



Anti-TDRP (aa 35-184) polyclonal antibody (DPABH-16367)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TDRP (testis development related protein) is a protein-coding gene.
Immunogen	antigen, corresponding to amino acids 35-184 of Human C8orf42 (Q86YL5).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [TDRP testis development related protein \[Homo sapiens \]](#)

Official Symbol	TDRP
Synonyms	TDRP; testis development related protein; Inm01; TDRP1; TDRP2; C8orf42; testis development-related protein;
Entrez Gene ID	157695
Protein Refseq	NP_001243042.1
UniProt ID	Q86YL5
Function	molecular_function;