



Anti-CCDC23 (aa 2-66) polyclonal antibody (DPABH-16373)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CCDC23 (coiled-coil domain containing 23) is a protein-coding gene. Diseases associated with CCDC23 include endotheliitis.
Immunogen	antigen: DPPARKEKTK VKESVSRVEK AKQKSAQQEL KQRQRAEIYA LNRVMTELEQ QQFDEFCKQM QPPGE, corresponding to amino acids 2-66 of Human CCDC23 (Q8N300).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	CCDC23 coiled-coil domain containing 24 [Homo sapiens]
Official Symbol	CCDC23
Synonyms	CCDC23; coiled-coil domain containing 23; coiled-coil domain-containing protein 23;
Entrez Gene ID	374969
Protein Refseq	NP_955374.1
UniProt ID	Q8N300