



# Anti-CSMD2 (aa 2462-2603) polyclonal antibody (DPABH-16422)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                            |                                                                                                                                                                    |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Antigen Description</b> | ESCO1 belongs to a conserved family of acetyltransferases involved in sister chromatid cohesion (Hou and Zou, 2005 [PubMed 15958495]).[supplied by OMIM, Mar 2008] |
| <b>Immunogen</b>           | antigen, corresponding to amino acids 2462-2603 of Human CSMD2                                                                                                     |
| <b>Isotype</b>             | IgG                                                                                                                                                                |
| <b>Source/Host</b>         | Rabbit                                                                                                                                                             |
| <b>Species Reactivity</b>  | Human                                                                                                                                                              |
| <b>Purification</b>        | Immunogen affinity purified                                                                                                                                        |
| <b>Conjugate</b>           | Unconjugated                                                                                                                                                       |
| <b>Applications</b>        | IHC-P                                                                                                                                                              |
| <b>Format</b>              | Liquid                                                                                                                                                             |
| <b>Size</b>                | 100 µl                                                                                                                                                             |
| <b>Buffer</b>              | pH: 7.20; Constituents: 59% PBS, 40% Glycerol                                                                                                                      |
| <b>Preservative</b>        | 0.02% Sodium Azide                                                                                                                                                 |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.                                                                              |

## GENE INFORMATION

**Gene Name** [CSMD2 CUB and Sushi multiple domains 3 \[ Homo sapiens \]](#)

|                        |                                                                                                                                                                                                       |
|------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Official Symbol</b> | CSMD2                                                                                                                                                                                                 |
| <b>Synonyms</b>        | CSMD2; CUB and Sushi multiple domains 2; dJ947L8.1; dJ1007G16.1; dJ1007G16.2; CUB and sushi domain-containing protein 2; CUB and Sushi (SCR repeat) domain; CUB and sushi multiple domains protein 2; |
| <b>Entrez Gene ID</b>  | <a href="#">114784</a>                                                                                                                                                                                |
| <b>Protein Refseq</b>  | <a href="#">NP_001268885.1</a>                                                                                                                                                                        |
| <b>UniProt ID</b>      | <a href="#">Q7Z408</a>                                                                                                                                                                                |