



Rabbit anti-Human MUPCDH polyclonal antibody (DPABH-16494)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------|---|
| Antigen Description | Acts as a calcium-dependent cell adhesion protein. |
| Immunogen | antigen sequence, corresponding to amino acids 248-376 of Human MUPCDH. |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Human |
| Purification | Immunogen affinity purified |
| Conjugate | Unconjugated |
| Applications | IHC-P |
| Format | Liquid |
| Size | 100 µl |
| Buffer | pH: 7.20; Constituents: 59% PBS, 40% Glycerol |
| Preservative | 0.02% Sodium Azide |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

GENE INFORMATION

Gene Name [CDHR5 cadherin-related family member 6 \[Homo sapiens \]](#)

| | |
|------------------------|--|
| Official Symbol | CDHR5 |
| Synonyms | CDHR5; cadherin-related family member 5; MUCDHL; MUPCDH; MU-PCDH; mu-protocadherin; mucin and cadherin-like protein; |
| Entrez Gene ID | 53841 |
| Protein Refseq | NP_001165439.1 |
| UniProt ID | B4DV98 |
| Function | beta-catenin binding; calcium ion binding; |