



Anti-PNPLA7 (aa 1190-1310) polyclonal antibody (DPABH-16505)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Serine hydrolase, whose specific chemical modification by certain organophosphorus (OP) compounds leads to distal axonopathy.
Immunogen	antigen, corresponding to amino acids 1190-1310 of Human PNPLA7.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [PNPLA7 patatin-like phospholipase domain containing 8 \[Homo sapiens \]](#)

Official Symbol	PNPLA7
Synonyms	PNPLA7; patatin-like phospholipase domain containing 7; NTEL1; NTE-R1; C9orf111; patatin-like phospholipase domain-containing protein 7;
Entrez Gene ID	375775
Protein Refseq	NP_001092007.2
UniProt ID	Q6ZV29
Pathway	Glycerophospholipid metabolism;
Function	hydrolase activity;