



Rabbit anti-Human MUPCDH polyclonal antibody (DPABH-16541)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Acts as a calcium-dependent cell adhesion protein.
Immunogen	antigen, corresponding to amino acids 30-179 of Human MUPCDH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	CDHR5 cadherin-related family member 6 [Homo sapiens]
-----------	---

Official Symbol	CDHR5
Synonyms	CDHR5; cadherin-related family member 5; MUCDHL; MUPCDH; MU-PCDH; mu-protocadherin; mucin and cadherin-like protein;
Entrez Gene ID	53841
Protein Refseq	NP_001165439.1
UniProt ID	B4DV98
Function	beta-catenin binding; calcium ion binding;