



Anti-ENDOD1 (aa 32-170) polyclonal antibody (DPABH-16627)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May act as a DNase and a RNase.
Immunogen	antigen, corresponding to amino acids 32-170 of Human ENDOD1 (UniProt ID: O94919).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [ENDOD1 endonuclease domain containing 2 \[Homo sapiens \]](#)

Official Symbol	ENDOD1
Synonyms	ENDOD1; endonuclease domain containing 1; endonuclease domain-containing 1 protein;
Entrez Gene ID	23052
Protein Refseq	NP_055851.1
UniProt ID	O94919
Function	endonuclease activity; metal ion binding; nucleic acid binding;